Suggested Electrode Placement Chart

The sticky electrode pads need to be placed in the correct position for maximum pain relief. It is preferable to obtain assistance from either your birthing partner or medical adviser when positioning the pads. Please read the following instructions carefully:

Bra line

Navel line

Panty line

of the skin.

Mar:4 (2): 50-2

1 1

2 2

Electrodes can be used several

times but the number of uses

PGCEA. Analgesia in Labour:

A review of the TENS method.

varies according to the condition

Reference: Cluett E. (1994) SRN RM ADM

Professional Care of Mother and Child.

Q Why does the stimulation not appear to be as strong

TENS 2 for a while?

setting.

after you have used the Elle

A Increase the intensity, you may

Q Why does the stimulation

A You may need to switch off

Q Why does the stimulation

even on a high setting?

O What do I do if there is no

stimulation, or weak

A Check electrode contact.

Call our customer careline

on +44 (0)20 8532 9595.

new battery?

Check leads.

Check batteries

stimulation intermittent

stimulation felt even with a

A The batteries may need

replacing.

electrode pads.

the unit and reposition the

not feel sufficiently strong

have become used to a lower

feel strong but ineffective?

1. Ensure unit remains switched off. Carefully peel the pads from the plastic backing using your forefinger and thumb. Do not pull on the lead wires.

2. Position the top two pads (connected to the left socket CH1) approximately 4cm either side of the spine and 4-5cm above the imaginary navel line as seen in the diagram.

3. Position the bottom two pads (connected to the right socket CH2) approximately 4cm either side of the spine and 4-5cm below the imaginary navel line as seen in the diagram.

Frequently Asked Questions

Q When shall I start using my Elle TENS 2?

A Use as early as you can after the onset of labour to allow time for your body's painkilling chemicals to rise.

Q Can I try the Elle TENS 2 before going into labour?

- A Yes, place 2 pads attached to one leadwire on your forearm to feel the different pulse sensations following the instructions on page 8.
- Q Can I combine the Elle TENS 2 with other medication?
- A Yes, you can still use analgesics such as gas and air (entinox) or pethidine
- O Can the Elle TENS 2 be used in hospital?
- A Yes. However, a midwife may ask you to switch it off temporarily.
- Q Can I use the Elle TENS 2 if I am considering a water birth?
- A Yes, but not in the water itself.

Using Your Elle TENS 2

a) You are now ready to start using vour Elle TENS 2. Switch the unit on by pressing the 🔘 on/off button. Ensure that the 'boost' button located at the top right of the unit is NOT pressed in. The Elle TENS 2 will always start up in BURST 1 mode.

b) Begin by pressing the 💽 CH1+ key. Each time you press this key, a small 'half-moon' shape will appear on the screen increasing in number and size as the intensity is increased. You will begin to feel a pulsating sensation which will burst in time with the small arrow flashing on the

screen. Ensure that you increase the intensity slowly. c) Repeat the previous step with the CH2+ key. You will now feel a

sensation in both sets of pads. d) BURST 1 (use between contractions): This is the first mode that you will use. It feels like a

bursting, tingling sensation i.e. the pulses are off and on periodically. This mode will help to promote those allimportant natural pain-killing chemicals known as 'endorphins' and 'encephalins'.

Troubleshooting Questions Useful Information

1. Continue stimulation for as long as necessary, unless you experience discomfort.

2. If you stimulate for a number of hours ensure the pads do not become dry. If this occurs, switch off the unit, disconnect the electrodes, and apply a small drop of

- 3. The pads can be used several times. Replace the pads when they stop sticking to the skin. Replacement pads are available from
- www.babycaretens.com. 4. Always ensure that the unit is switched off before you
- remove all the leads and pads. Do not pull the wires, ALWAYS remove the jack plug from its socket by holding the plug between the forefinger and thumb and pulling gently.
- 5. The Elle TENS 2 can also be used as a standard TENS unit for relief of numerous conditions, including after pains and post operative pains It is worth bearing this in mind if you have a caesarean section. Please seek medical

advice first.

e) BOOST 1 (use during

continuous sensation.

contractions): When experiencing a contraction, simply press the 'boost' button located at the top right of the unit. The Elle TENS 2 will then switch automatically into the BOOST 1 mode for that extra surge of power required to combat the pain during contractions. This will feel like a

f) Once the contraction has passed,

simply release the 'boost' button. The

Elle TENS 2 will return to the original

BURST 1 mode. Continue this cycle

during your early stages of labour.

Using the Contraction Timer

automatically trigger the Contraction

Timer (CT). This can be seen in the

bottom right of the display. This

measures the duration of your

Button is released, the unit will

trigger the Between Contraction

Timer (BCT) which measures the

be seen in the bottom left of the

display.

Diagrams

KEY:

time between contractions. This can

Post-Natal Pad Placement

contractions in seconds. Once the

contraction has passed and the Boost

The first time you press the Boost

Button, the Elle TENS 2 will

history To recall a list of your last 5 'Contraction Times' (CT) and

'Between Contraction Times' (BCT) press the 🕒 🔗 CH1- and CH2buttons simultaneously. The last recorded set of BCT and CT times are shown as 1.

Recalling your contraction

Use 🕜 CH2+ to go back to previous sets (eg, 2, 3, 4, 5) or forward again by pressing SCH2key. Recalling the contraction history will not interrupt the functionality of the unit. N.B. To avoid losing contraction history, do NOT switch unit off. To return back to the original mode screen, (which will have continued counting in the background) press the 🕒 🖉 CH1- and CH2- buttons together again.

The Mode Button

There are two mode ranges available on the Elle TENS 2, each with a Burst and Boost function. Mode 1 (BURST 1 and BOOST 1), which has been described above, is primarily used for the earlier stages

of labour. Mode 2 (BURST 2 and BOOST 2) usually used when your labour

Elle TENS 2 Technical Data

Programme Display Output Frequency / Duration	Burst Burst hund per second, 8 pulse record a pulse record. 8 pulser per burst*	en BURST Burst Burst Publie Frequency (32H3) Pulse Width (150µs) 2 bursts per second, 8 pulses per burst*	Boost Boost 1 Pulse Frequency (80H4) Pulse Width (150µs)		Boost 2 Boost 2 ruse requency (work frage what i Lage	d Output Channel Dual Channel		Output Waveform Symmetrical bi-phasic rectangular	Output Voltage 0 - 55v zero to peak adjustable in 15 steps at 5000 load	Output Intensity 0 - 110mA adjustable in 15 steps at 5000 load Intensity level will dop back to see after programme changes	 Pulse width adjustment does not charge when switching between Bust and Boost modes. Delasting width is Styles. There are 5 adjustable pulse widths selectable in all modes. Storm and Storm an
	Equipment type Dual Channel DigitalTENS Dual Channel DigitalTENS This unit is fully compliant with FEC Directive 93/42/FEC, This unit is fully compliant with FEC Directive 93/42/FEC, tasking a internation proved equipment type BF. The Interded for compliance equipment type BF.			75 grams	CE Touch-proof jacks	Digital LCD - 8 seven segment digits and around 95 segments	Will switch off after 5 minutes if not used	22v ± 02v	2 × 1.5v AA Battery (Type LR6)	Temperature 0-40°C Air humidity 30-75% Air pressure 700 hPa-1060 hPa	Temperature 040°C Air humidity 30-90% Air pressure 700 hPa- 1060 hPa
Equipment type	Dual Channel Digital TENS Classification This unit is fully compliant classified as internally pow	It is intended for continuous operation. *Default pulse width after mode changed (pube width no and 8005T mode) 5 selectable pulse widths in all modes.	Unit Size	Unit Weight	Output Jacks	Display	Automatic Power Off	Battery Low Detect	Battery	Operating Environment	Storage and Transport Environment

advances and your contractions become more frequent.

Press the 🙆 Mode Button to switch between mode ranges. You will see from the User Display, that the number on the centre right of the display will change from 1 to 2 (and vice-versa). N.B. Both modes are interchangeable and can be used according to your preference, not necessarily for the stage of labour you are in.

The Opti-Max Key

This provides extra versatility to your Elle TENS 2. Increasing or decreasing the 🙆 Opti-Max Key will allow you to either heighten or lessen the level of pulse power transmitted to the electrode pads. Pressing the Opti-Max Key either + or - will increase or decrease Opti-Max sensation. This is extremely useful for the final stages of labour when you need an extra surge of power. There are 5 Opti-Max levels: 50, 100, 150, 200 and 250µs (microseconds). The unit will always start up at Opti-Max 150 µs. which is exactly half way between the highest and lowest setting. N.B. The Elle TENS 2 starts up in BURST 1. The BCT only appears after one cycle of pressing the Boost button on and off.





Guarantee

Your Elle TENS 2 machine is guaranteed for a period of 2 years against manufacturer's defects. Please record purchase details in the space provided below for your record purposes. The guarantee does not include leads, electrodes or battery. Applies to purchases only.

Date of purchase:

Invoice number:




Body Clock Health Care Ltd 108 George Lane, South Woodford, London E18 1AD

Tel (+44) (0)20 8532 9595 Fax (+44) (0)20 8532 9551 www.babycaretens.com e-mail customercare@babycaretens.com

> Babycare TENS is part of the Body Clock group of companies

 \mathbf{E} 0120 Jun 2018 (2)

recycle This manual is recyclable



= first pair of electrodes







Outer and inner aspect of shoulder

water to them.

Elle TENS 2 **Quick Reference Guide**


If you are in labour, please use the quick reference guide set out below. It is however recommended that you read the section entitled "Setting Up Your Elle TENS 2" found on page 5 to obtain the best performance from your unit.

Before using the Elle TENS 2, please read the contraindications and warnings on pages 1 and 2.

- 1. Ensure the Elle TENS 2 unit is switched off.
- 2. Insert 2 x AA batteries in the compartment at the back of the unit (see figs 1 and 2, page 6).
- 3. Insert the leadwire plugs into the top of the unit (see fig 3, page 6).
- 4. Insert the pins at the other end of the leadwires into the sockets attached to the electrode pads. Do not remove the electrodes from their backing at this stage (see fig 4, page 6).

Controls For The Elle TENS 2

instructions to obtain the most out of the unit.



5. Position the electrode pads on Contents

Elle TENS 2 Quick Reference Guide

General Precautions....

Useful Diagrams.....

on Unit

Λ

Using Your Elle TENS 2......

How Will The Elle TENS 2 Help Your Pain?.... 1

Introduction To The Elle TENS 2 1

Precautions and Contraindications......1

Controls For The Elle TENS 2...... 3

Display Symbols and their Meanings...... 4

Setting Up Your Elle TENS 2 5

Frequently Asked Questions 10

Post Natal Pad Placement Diagrams......11

Equipment providing a

protection against electric

allowable leakage currents

having an F type (floating)

Warning - refer to page 2 of

shock particularly regarding

particular degree of

Explanation of Symbols

applied part.

these instructions

6

8

the body as indicated on page 7. A pair of electrodes must be connected to each lead wire in sensation. The electrode pads must be peeled away from their backing and placed directly on your body for you to feel any

- 6. Make sure that the 'boost'
- BURST 1 mode (see page 3).
- 8. Now push the 💽 (CH1 +) key until you feel a pulsating sensation.
- 9. Repeat the same with the (CH2+) key.
- - the contraction has passed and continue in the BURST 1 mode as before

The Elle TENS 2 is very simple to use. However, it is essential to read these

order for you to feel any sensation.

- button at the side of the unit is NOT pressed in. This ensures that the unit will start in the
- 7. Press the (0) (on/off) to turn the unit on.

- 10. At the onset of a contraction press the 🗲 'boost' button for the additional surge of power needed to combat the pain.
- 11. Release the 'boost' button once

Display Symbols And Their Meanings

BURST

BURST is used at the onset of labour and in between contractions. There are two Burst modes available indicated by a 1 or 2. You can switch between them by pressing the Mode Button (see page 9).

►BOOST

The Boost Button at the top right of the unit activates the BOOST mode This function is used during contractions. There are two Boost modes available indicated by a 1 or 2. You can switch between them by pressing the Mode Button.

BETWEEN CONTRACTION TIMER Starts by pressing the Boost Button, this automatically measures the time interval between your contractions.

▶BURST 1 ▶BURST 2

MODE 1 AND 2 There are two mode ranges available each with a Burst and Boost function (see page 9).



IP22 The first number 2: Protected against access to hazardous parts with a finger, and the jointed test finger of 12 mm\$\Phi, 80 mm length, shall have adequate clearance from hazardous parts, and protected against solid foreign objects of 12.5 mm ϕ and greater.

The second number 2: Protected against vertically falling water drops when enclosure tilted up to 15° Vertically falling drops shall have no harmful effects when the enclosure is tilted at any angle up to 15° on either side of the vertical.



This symbol indicates the serial number of the device and includes the year of manufacture. The serial number can be found in the battery compartment of this

This unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENT.

SN

device.

LOW BATTERY INDICATOR

This sign will start flashing if the

batteries run low. If this is the case,

change both batteries. Please refer

This shows the output from 0 to 15

Starts by pressing the Boost Button,

this automatically records the

duration of your contractions.

Indicates the Opti-Max output.

Refer to page 9 for more details.

This symbol will appear if any of

your pads are not attached securely

PAD CONTACT WARNING

Each segment represents approx

to page 5 and 6 for instructions.

for both channel 1 and 2.

7% of the total output.

CONTRACTION TIMER

Max ull 250us

OPTI-MAX

-12---


to the skin

(safety cut-out)

POWER BAR

GDA

BCD ALA



How Will The Elle TENS 2 Help Your Pain?

TENS (Transcutaneous Electrical Nerve Stimulation) works by passing mild electrical impulses through the skin, via electrode pads, into the nerve fibres which lie below. The TENS impulses help your body produce its own pain killing chemicals, such as 'endorphins' and 'encephalins'. The level of pain relief obtained varies from person to person. Some find that the Elle TENS 2 provides all the pain relief required during the whole of labour whereas others find that other analgesics are also

Introduction To The Elle TENS 2

labour

needed during the later stages of

The Elle TENS 2[™] is an efficient method of drug-free pain relief. This powerful digital TENS unit is smooth, sophisticated and simple to use. It has been cleverly crafted by Babycare TENS specifically for home use for mums-to-be, for the management of pain for labour and beyond.

Setting Up Your Elle TENS 2

Check the following contents:

- 1 x Elle TENS 2 Unit
- 1 x Pack of 4 Self Adhesive Electrodes (40mm x 100mm)
- 2 x Leadwires and 1 Spare
- 2 x AA Batteries and 2 Spare 1 x Carrying Pouch
- 1 x Easy-Release Neck Cord
 - 1 x Instruction Manual
- 1 x EMC documentation

The Assembly Stage

- Make sure the unit is swiched off
- a) Press and rotate the belt clip at the back of the Elle TENS 2 in either direction to reveal the battery compartment.
- b) Release the battery compartment cover and insert 2 x AA batteries ensuring the positive (+) and negative (-) terminals are correctly positioned as marked in the battery compartment. Replace the cover.



- 1. TENS must NOT be used before the 37th week of pregnancy.
- 2. If you are pregnant, do NOT place the electrodes over your abdomen.
- 3 Do NOT use if you have a pacemaker or any other implanted electrical device.
- Check with your healthcare 4 professional before using if you suffer from a heart condition, epilepsy, undiagnosed pain, have any metal implants, or any doubts whatsoever.
- 5. Do NOT use if you experience an allergic reaction to the electrodes
- 6 Never use TENS to mask undiagnosed pain since this could require urgent treatment.

Take two of the purple

TENS 2 unit.

electrodes

black pin. It makes no

they are attached to.

e) If you wish, attach the easy-

release neck cord which

TENS 2 from your neck.

Thread the loop of the neck cord

through the space at the top of

neck cord through the loop and

the unit. Slip the other end of the

Using the Neck Cord for

mobility during labour

pull into position.

enables you to hang the Elle

difference which electrode

leadwires and insert each of

the black jack plugs into the


sockets at the top of the Elle

c)

General Precautions

- 1. Do NOT use this unit without first reading these instructions.
- 2. Do NOT immerse the Elle TENS 2 in any liquid. Do not use in the bath or the shower.
- 3. Do NOT place it close to any source of excessive heat or operate it in the presence of flammable gas.
- 4. Do NOT drop this unit onto a hard surface.
- 5 Do NOT attempt to dismantle the Elle TENS 2.
- 6. Only use specified batteries and electrodes.
- 7. If damaged, do not use. Return to supplier.
- Remove the batteries when 8. not in use
- 9. Do NOT use while driving or operating potentially dangerous machinery or while using a microwave.
- 10. Keep out of the reach of children.

Useful Diagrams



or near the eyes, in the mouth, over the front or sides of the neck or across the head or the heart.

11. Do NOT place electrodes on

- 12. Do NOT place the electrodes on varicose veins or recent scarring or on any area of broken, inflamed, infected or numb skin. Electrodes should generally only be applied to skin with normal sensation unless under medical supervision.
- 13. Do NOT use in the presence of tuberculosis, malignant tumours, very high or very low blood pressure, high fever or acute inflammatory disease unless under medical supervision.
- 14. A mild shock may be caused if one electrode is accidentally removed in the event of this happening the unit should be turned off immediately.

NB. The unit should be wiped clean periodically using medically approved anti-bacterial wipes or a very slightly dampened cloth. A mild soap may be used but DO NOT apply solvents.

d) Now insert the pins (at the other end of the leadwires) into the sockets attached to the sticky Do not remove the electrodes from the backing at this stage. NB. Each leadwire has a red and

Sticky Electrodes